

Digitales Catering-Thermometer Typ T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Catering-Thermometer

Hochgenaues Kochthermometer

Thermoelement Typ T

HACCP-konform,

-100 bis + 280C

Hold-Funktion,

IP67 wasserdichtes Thermometer

Duales LCD-Display

Einfache Bedienung: Zwei-Tasten-Bedienung

Duales Display

Automatische Abschaltfunktion.

Funktionen:

Selbstkalibrierend

Anzeige für schwache Batterie

Lieferung mit 1 9V PP3 Batterie

Messbereich -100 bis +280C

Hohe Genauigkeit +/-0,1 % des Messwerts

+/-0,2 C, 0,1 C Auflösung

Der CA2005 bietet eine Hold-Funktion und wenn diese Taste gedrückt wird, wird die aktuelle Temperatur gehalten und an die obere rechte Ecke des Displays übertragen. Das Hauptdisplay zeigt weiterhin die Temperatur in Echtzeit an.

Leistungsbeschreibung

Specifications

Product Code	XE-4090-001
General Description	High accuracy chef thermometer
Type	IP67 waterproof thermometer



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Accuracy	+/-0.1% of Reading
Resolution	+/-0.2°C, 0.1°C Resolution
Input / Ranges	-100 to + 280°C
Power Supply	Supplied with 1 9V PP3 battery
Sensor Type	T