



**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## PTFE-Single-Shot-Thermoelementkabel Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**PTFE-Single-Shot-Thermoelementkabel Typ K**

Einzigartige Single-Shot-Kabelkonstruktion Beständig gegen Chemikalien, Öle und Flüssigkeiten. Geeignet für Anwendungen in der Luft- und Raumfahrt und wird häufig in Autoklaven verwendet.

Spezifikationen

## Specifications

<b>General Description</b>	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.
<b>Cable Length</b>	10, 25, 50, 100M
<b>Cable Type</b>	PTFE
<b>Core / Strands</b>	1/0.376mm
<b>Colour</b>	Green
<b>Max. Temperature</b>	+250°C
<b>Standards Met</b>	IEC-584-3
<b>Thermocouple Type</b>	K