



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PTFE-Single-Shot-Thermoelementkabel Typ T



50M



10m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PTFE-Single-Shot-Thermoelementkabel Typ T

Einzigartige Single-Shot-Kabelkonstruktion Beständig gegen Chemikalien, Öle und Flüssigkeiten. Geeignet für Anwendungen in der Luft- und Raumfahrt und wird häufig in Autoklaven verwendet.

Spezifikationen

Specifications

General Description	Unique single shot cable construction Resistant to chemicals, oils and fluids. Suitable for aerospace applications and is widely used in autoclaves.
Cable Length	10, 25, 50, 100M
Cable Type	PTFE
Core / Strands	1/0.376mm
Colour	Brown
Max. Temperature	+250°C
Standards Met	IEC-584-3, Class 1
Thermocouple Type	T