



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Pt100 Klasse B Sensoren mit teflonisolierter Leitung - Typ PRT - 4,0 mm Durchmesser - 125 mm Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

Product Code	XE-3680-001
Sensor Type	Pt100 element in a 3.0mm or 4.0mm diameter 316 stainless steel probe
Element	Pt100 element to IEC 751 Class B
Probe Diameter	4.0mm
Probe Length	125mm
Sheath Material	316 stainless steel probe
Cable Length	1M
Cable Type	Teflon® insulated twisted lead
Core / Strands	7/0.16mm
Cable Termination	Bare tails
Max. Temperature	200°C
Min. Temperature	-50°C
Standards Met	IEC 751
Accuracy	Class B