

## Pt1000 Klasse-B-Sensoren mit teflonisoliertem Kabel - Typ PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Insulated element in a 4.0mm diameter 316 stainless steel probe
<b>Sensor Type</b>	Pt1000 (1000 Ohms @ 0°C) to IEC 751, Class B
<b>Element</b>	Pt1000
<b>Probe Diameter</b>	4mm
<b>Sheath Material</b>	316 stainless steel probe
<b>Cable Length</b>	1m
<b>Cable Type</b>	Teflon® insulated twisted leads
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-50°C
<b>Number of Wires</b>	2