



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Pt1000 Klasse B Sensoren mit teflonisolierter Leitung - Typ PRT - 4,0 mm Durchmesser - 125 mm Länge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Leistungsbeschreibung

**Specifications**

<b>Product Code</b>	XE-3684-001
<b>General Description</b>	Insulated element in a 4.0mm diameter 316 stainless steel probe
<b>Sensor Type</b>	Pt1000 (1000 Ohms @ 0°C) to IEC 751, Class B
<b>Element</b>	Pt1000
<b>Probe Diameter</b>	4.0mm
<b>Probe Length</b>	125mm
<b>Sheath Material</b>	316 stainless steel probe
<b>Cable Length</b>	1M
<b>Cable Type</b>	Teflon® insulated twisted leads
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	200°C
<b>Min. Temperature</b>	-50°C
<b>Number of Wires</b>	2
<b>Series</b>	Pt1000 Sensors