


Typ J Glasfaser isoliert Flachpaar Thermoelement Kabel / Draht (BS)

Type J, 1/0.508mm
Fiberglass Flat pair



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typ J Glasfaser isoliert Flachpaar Thermoelement Kabel / Draht (BS)

Bild ist nur für Farbcodierung summieren

Der Preis ist pro Meter - bitte fügen Sie die erforderliche Menge für den Korb. Nur noch 1.227 Meter auf Lager

Glasfaser isolierte Flachpaarkonstruktion, Leiter flach gelegt, Glasfaser mit Glasfaser insgesamt isoliert, Silikon lackiert. Gute Temperaturbeständigkeit, aber nicht geeignet, wo Flüssigkeiten vorhanden sind.

Spezifikationen
Specifications

Product Code	WJ-033-D
General Description	Good temperature range but will not prevent ingress of fluids. Fairly flexible but does not provide good mechanical protection.
Cable Length	Priced per metre
Cable Type	Flat pair cable
Core / Strands	1/0.508mm
Outer Colour	Black
+ Colour	Yellow
- Colour	Blue
Insulation Material	Glassfibre insulated
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
+ Conductor	Iron
- Conductor	Constantan
Standards Met	BS



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple Type	J
Tolerance	Class 2