


Thermoelementkabel / Draht Typ K PVC isoliert flach paar (BS)

Type K, 13/0.2mm
PVC Flat Pair



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelementkabel / Draht Typ K PVC isoliert flach paar (BS)

Bild ist nur für Farbcodierung summieren

Der Preis ist pro Meter - bitte fügen Sie die erforderliche Menge für den Korb. Nur noch 613 Meter auf Lager

Das Kompensationskabel ähnelt der thermoelektrischen Eigenschaft des relevanten Thermoelements über einen begrenzten Umgebungstemperaturbereich, typischerweise 0° bis 80°C. Da diese Kabel aus kostengünstigen Materialien hergestellt werden, können bei Anlagenanlagen Kosteneinsparungen erzielt werden, verglichen mit dem laufenden echten Thermoelement-Verlängerungskabel.

Spezifikationen
Specifications

Product Code	WK-010-D
General Description	Good general purpose insulation for 'light' environments. Waterproof and very flexible.
Cable Length	Priced per metre
Cable Type	Flat pair construction
Core / Strands	13/0.2mm
Outer Colour	Red
+ Colour	Brown
- Colour	Blue
Insulation Material	PVC Insulated
Max. Temperature	+105°C
Min. Temperature	-10°C
+ Conductor	Nickel Chromium
- Conductor	Nickel Aluminium



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standards Met	BS
Thermocouple Type	K
Tolerance	Class 2