


# Miniatur-Thermoelement-Stecker FMTC-CU-M Typ Cu



**Miniature in-line Plug**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Thermoelement-Anschlussstecker FMTC-CU-M Typ CU**

Mit flachen Stiften, die dem Thermoelementtyp entsprechen und sich ideal für den allgemeinen Gebrauch eignen, sind alle Kontakte polarisiert, um eine korrekte Verbindung zu gewährleisten. Diese Steckverbinder nehmen Thermoelementkabel mit einem Durchmesser von bis zu 4 mm auf.

**Leistungsbeschreibung**
**Specifications**

<b>Product Code</b>	XE-1370-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Flat Pins
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS