

## Miniatur-Thermoelement-Steckverbinder Inline-Buchse IM-J-FS Typ J IEC



**Miniature in-line Socket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Thermoelement-Steckverbinder Inline-Buchse IM-J-FS Typ J IEC**

Ideal fr den allgemeinen Gebrauch, alle Kontakte sind polarisiert, um eine korrekte Verbindung zu gewhrleisten.  
Diese Steckverbinder nehmen Thermoelementkabel mit einem Durchmesser von bis zu 4 mm auf.

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1162-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	J IEC
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Black
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Iron
<b>-Negative Contact</b>	Copper Nickel
<b>Standards Met</b>	BS EN 50212:1997. RoHS