





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Miniatur-Thermoelement-Steckverbinder Inline-Buchse IM-N-FS Typ N IEC



**Miniature in-line Socket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Thermoelement-Steckverbinder Inline-Buchse IM-N-FS Typ N IEC**

Ideal for general use, all contacts are polarised, to ensure correct connection.  
These connectors take thermocouple cable with a diameter of up to 4 mm.

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1252-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	N IEC
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Pink
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel / Chromium / Silicon
<b>-Negative Contact</b>	Nickel Silicon
<b>Standards Met</b>	BS EN 50212:1997. RoHS