


Verstellbares Bajonett, berflochtenes Glasfaserkabel aus Edelstahl - Typ PRT - Lnge 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Verstellbares Bajonett mit 200 mm Feder mit verstellbarer Befestigungskappe (11,5 mm Innendurchmesser Einzelschlitz).

Leistungsbeschreibung

Specifications

Product Code	XE-5655-001
General Description	PRT with adjustable insertion depth
Sensor Type	Pt100 4-wire
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel Overbraid
Core / Strands	4 core 7/0.2mm
Cable Termination	Bare tails
Max. Temperature	350°C
Min. Temperature	-50°C
Standards Met	IEC-60751 Class B
Bayonet cap	11.5mm I.D. single slot
Spring	6.0mm diameter x 200mm long