


Einstellbares Bajonett Pt100, Glasfaser Edelstahl Overbraided Kabel - Typ RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Einstellbares Bajonett PRT.

Einstellbares Bajonett (Pt100) PRT / RTD mit 200mm Feder mit verstellbarer Befestigungskappe (11,5 mm I.D. Einzelsteckplatz)

Erhältlich in 1M oder 2M Längen.

Specifications

General Description	PRT with adjustable insertion depth
Sensor Type	Pt100 4 wire
Cable Length	1 or 2 metres
Cable Type	Glassfibre insulated lead with stainless steel Overbraid
Core / Strands	4 core 7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Min. Temperature	-50°C
Standards Met	IEC-60751 Class B
Bayonet cap	11.5mm I.D. single slot
Spring	6.0mm diameter x 200mm long