

Miniatur-Schnellkupplungs-Thermoelement-Buchse FMTC-CU-FQ Typ Cu



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Schnellkupplungs-Thermoelement-Buchse FMTC-CU-FQ Typ CU

Einfache Jab-In-Verbindung, einfach den Draht eindrcken und die Schraube festziehen. Die Kontakte sind polarisiert, um einen falschen Anschluss zu verhindern.

Leistungsbeschreibung**Specifications**

Product Code	XE-1382-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	White
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS