

Miniatur-Rechteckiger Thermoelement-Steckverbinder Blendenbuchse IM-N-FF Typ N IEC



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Rechteckiger Thermoelement-Steckverbinder Blendenbuchse IM-N-FF Typ N IEC

Nylon-Clip zur schnellen und einfachen Befestigung an der Platte. Alle Blendenmontagebuchsen sind mittels ungleich groer Stifte polarisiert

Leistungsbeschreibung**Specifications**

Product Code	XE-1257-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	N IEC
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Pink
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS