



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Pt100 Duplex 2 x 3 Leiter Klasse A Widerstandsthermometer, KNE-Kopf - 150 mm Sondenlänge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

Product Code	XE-3656-001
Sensor Type	2 x Pt100 (100 Ohms @ 0°C) to IEC 751, Class A
Probe Diameter	6mm
Probe Length	150mm
Sheath Material	316 stainless steel
Head Type	KNE
Termination Type	6 way connection block for KNE head
Max. Temperature	+250°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class A
Number of Wires	2 x 3 Wire
Series	Pt100 KNE