

## **Pt100 Duplex 2 x 3 Leiter Klasse A Widerstandsthermometer, KNE-Kopf - 200 mm Sondenlnge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Leistungsbeschreibung

**Specifications**

<b>Product Code</b>	XE-3657-001
<b>Sensor Type</b>	2 x Pt100 (100 Ohms @ 0°C) to IEC 751, Class A
<b>Probe Diameter</b>	6mm
<b>Probe Length</b>	200mm
<b>Sheath Material</b>	316 stainless steel
<b>Head Type</b>	KNE
<b>Termination Type</b>	6 way connection block for KNE head
<b>Max. Temperature</b>	+250°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP67
<b>Accuracy</b>	Class A
<b>Number of Wires</b>	2 x 3 Wire
<b>Series</b>	Pt100 KNE