



Miniatur Thermoelement Stecker & Socket AM-B-M+F Typ B ANSI



Miniatur in-line Plug
(dimensions in mm)



Miniatur in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-1017-001
General Description	Ideal for general purpose use. Contacts are polarised to ensure correct connection
Thermocouple Type	B ANSI
Connector Size	Miniature
Connector Type	Plug & Socket
Colour	Grey
Max. Temperature	220°C
Pins	Flat Pins
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS