


Miniatur-Schnelldraht-Thermoelement-Steckverbinder Stecker & Buchse AM-K-MQ+FQ Typ K ANSI



Miniature in-line Plug - Quick Wire
(dimensions in mm)



Miniature in-line Socket - Quick Wire
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Steckverbinder Schnelldrahtstecker & Buchse AM-K-MQ+FQ Typ K ANSI

Miniatur-Steckverbinder mit flachen Steckstiften. Einfacher Jab-in-Anschluss fr eine schnelle Terminierung. Einfach Draht einstecken & Schraube festziehen.

Leistungsbeschreibung**Specifications**

Product Code	XE-1020-001
General Description	Easy jab-in connection for rapid termination
Thermocouple Type	K ANSI
Connector Size	Miniature
Connector Type	Quick Wire Plug & Socket
Colour	Yellow
Max. Temperature	220°C
Pins	Flat Pins
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector