


Standard-Thermoelement-Steckverbinder Stecker & Buchse IS-K-M+F Typ K IEC



Standard in-line Plug
(dimensions in mm)



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelement-Steckverbinder Stecker & Buchse IS-K-M+F Typ K IEC

Die Kontakte werden aus echten Thermoelement-Legierungen hergestellt. Geeignet für industrielle Anwendungen

Leistungsbeschreibung**Specifications**

Product Code	XE-1042-001
General Description	Contacts are polarised to prevent incorrect connection
Thermocouple Type	K IEC
Connector Size	Standard
Connector Type	Plug & Socket
Colour	Green
Max. Temperature	220°C
Pins	Solid Round Pins
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Standard Connector