

Mineralisiertes Thermoelement mit 0,25 mm Durchmesser und ANSI-Miniaturstecker - Typ K



IMAGE
COMING
SOON

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bitte verwenden Sie unsere [kundenspezifischen Sensorformulare](#) , um Ihre Anforderungen für verschiedene MI-Längen zu ermitteln, da alle mit einem IEC-einen schnellen und einfachen Anschluss abgeschlossen sind.

Müssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rückführbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rückführbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

Bitte beachten Sie, dass Sie bezüglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufügen.

Leistungsbeschreibung

Specifications

| | |
|------------------------------|--|
| General Description | Highly flexible, sheath can be bent/formed to suit many applications and processes |
| Sensor Type | Type K ANSI MB STYLE probes 0.25mm |
| Thermocouple Junction | Single element, junction insulated or grounded |
| Probe Length | 100mm / 150mm / 250mm |
| Sheath Material | Stainless steel 316 sheath |
| Termination Type | Fitted Miniature Plug Type ANSI |
| Range | -40 to +260°C |
| Important Information | A 1%+ Length Tolerance applies |