



Leiterplattenmontage Thermoelement-Steckverbinderbuchse IM-K-PCB Typ K IEC



Miniature PCB Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leiterplattenmontage Thermoelement-Steckverbinderbuchse IM-K-PCB Typ K IEC

Speziell fr Leiterplatten, massive Pins fr die direkte Platinenmontage und Bereitstellung fr CJC-Sensor

Leistungsbeschreibung**Specifications**

Product Code	XE-1109-001
General Description	Versatile, high quality circuit board mounting miniature socket
Thermocouple Type	K IEC
Connector Type	PCB Socket
Connector Size	Miniature
Colour	Green
Max. Temperature	220°C
Pins	Brass
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS