



# Leiterplattenmontage Thermoelement-Steckerbuchse IM-T-PCB Typ T IEC



**Miniature PCB Socket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leiterplattenmontage Thermoelement-Steckerbuchse IM-T-PCB Typ T IEC**

Speziell fr Leiterplatten, massive Pins fr die direkte Platinenmontage und Bereitstellung fr CJC-Sensor

**Leistungsbeschreibung****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1199-001   |
| <b>General Description</b> | Versatile, high quality circuit board mounting miniature socket |
| <b>Thermocouple Type</b>   | T IEC   |
| <b>Connector Type</b>      | PCB socket  |
| <b>Connector Size</b>      | Miniature   |
| <b>Colour</b>              | Brown   |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper Nickel   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |