

## Lascar EL-P-VACX High accuracy Glycol thermistor probe, uncalibrated



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

High accuracy, uncalibrated 3 metre long thermistor probe in a sealed Glycol-filled bottle for vaccine monitoring. For use with the EL-WiFi-VACX data logger

[\(XE-4058-001\)](#)

### Specifications

## Specifications

<b>Product Code</b>	XE-4082-001
<b>Model</b>	EL-P-VACX
<b>Type</b>	High accuracy Glycol thermistor probe, uncalibrated
<b>Accuracy</b>	+/-0.1°C / +/-0.2°F
<b>Cable Length</b>	3 metres
<b>Measuring range</b>	-40 to +60°C / -40 to + 140°F