

## 0,75 mm Durchmesser mit Miniatur-IEC-Stecker isoliert oder geerdet



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Typ K Mineralisolierter Edelstahl 316, Ausgestattet mit Ministecker

Wenn Sie grere Mengen bentigen, senden Sie bitte eine E-Mail [sales@labfacility.com](mailto:sales@labfacility.com) oder rufen Sie unsere Vertriebs hotline 01243 871280 an.

### **Mssen Sie Ihren Temperatursensor kalibrieren?**

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

**Bitte beachten Sie, dass Sie bezglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufgen. .**

### **Leistungsbeschreibung**

## **Specifications**

<b>General Description</b>	Type K 316 stainless steel, Fitted with miniplug
<b>Sensor Type</b>	Mineral Insulated
<b>Thermocouple Type</b>	K
<b>Element/hot junction</b>	Single element, junction insulated or grounded from sheath
<b>Probe Diameter</b>	0.75mm
<b>Probe Length</b>	50, 100, 150, 250mm
<b>Sheath Material</b>	316 stainless steel
<b>Termination Type</b>	Thermocouple IEC Mini plug (Max Temp 200°C)
<b>Max. Temperature</b>	+750°C
<b>Min. Temperature</b>	-100°C
<b>Tolerance</b>	Class 1
<b>Standards Met</b>	In accordance with IEC 584