


Digitales Haushaltsthermometer Hygrometer mit rückseitigem Magneten



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TypT-1/0.2

Specifications

Product Code	XE-4352-001
Model	LDTH-8
General Description	Digital Household Thermometer Hygrometer with rear magnet
Type	Timer-Thermometer 2 in 1
Accuracy	+/- 3%
Readings	Timer range: 99 Min 59 Sec
Resolution	0.1°C(0.2°F)
High Alarm Level	High temperature alarm function
Battery	CR2430
Cable Length	Wire: 120cm length
Length	10m
Colour	White
Dimensions	118mm x 60mm x 7.5mm
Functions	°C/°F
Function #1	Continuous: 33°F- 140°F
Display Type	LCD
Operating Humidity	25%-98% RH
Temperature Range	-50-70°C (-58-158°F)
Probe Material	304 stainless steel
Weight	86g
Sensor Type	Sensirion SHT30



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Standards Met	CE, FCC, ROHS
Construction	Eco-friendly ABS plastic + Stainless steel front shell
Display	FSLCD
Important Information	Support stand or magnetic mounting options