

## Digitales Thermometer mit klappbarem Fhler LDT802



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Product Code</b>	XE-4351-001
<b>Model</b>	LDT802
<b>General Description</b>	Digital folding probe thermometer
<b>Accuracy</b>	+/-1°C ( -2°F ) at -20°C to150°C (-4°F to 392°F)
<b>Readings</b>	Min max temperature record
<b>Resolution</b>	0.1°C/0.2°F
<b>Response</b>	Super fast response time of 3-4 seconds
<b>Battery</b>	3V button battery (included)
<b>Display Type</b>	Bright backlit
<b>Automatic Shutdown</b>	Power saving: 10 minutes auto-off
<b>Temperature Range</b>	-50°C to 300°C (-58°F to 572°F)
<b>Probe Diameter</b>	Foldable probe Diameter: 3.5mm
<b>Probe Material</b>	304 stainless steel
<b>Standards Met</b>	IP67 Waterproof design
<b>Construction</b>	Eco-friendly ABS plastic
<b>Important Information</b>	Inside magnet on backside