

Digitales Thermometer mit klappbarem Fhler LDT802



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Product Code	XE-4351-001
Model	LDT802
General Description	Digital folding probe thermometer
Accuracy	+/-1°C (-2°F) at -20°C to150°C (-4°F to 392°F)
Readings	Min max temperature record
Resolution	0.1°C/0.2°F
Response	Super fast response time of 3-4 seconds
Battery	3V button battery (included)
Display Type	Bright backlit
Automatic Shutdown	Power saving: 10 minutes auto-off
Temperature Range	-50°C to 300°C (-58°F to 572°F)
Probe Diameter	Foldable probe Diameter: 3.5mm
Probe Material	304 stainless steel
Standards Met	IP67 Waterproof design
Construction	Eco-friendly ABS plastic
Important Information	Inside magnet on backside