



Standard-Thermoelement-Stecker IS-N-M Typ N IEC



Standard in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelement-Stecker IS-N-M Typ N IEC

Die Kontakte werden aus echten Thermoelement-Legierungen hergestellt. Geeignet für industrielle Anwendungen

Leistungsbeschreibung**Specifications**

Product Code	XE-1550-001
General Description	Contacts are polarised to prevent incorrect connection
Thermocouple Type	N IEC
Connector Type	Plug
Connector Size	Standard
Colour	Pink
Max. Temperature	220°C
Pins	Solid Round Pins
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS