



**IEC Typ K 1,9 kg
Zugknopf-Magnet-Thermoelement, PFA-isoliertes
Kabel mit Edelstahl-bergeflecht mit Miniatur- oder
Standardstecker**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	1.9kg Pull Button Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	Type K
Magnet	19.1mm Diameter / 12.7mm height. Max Pull: 1.9kg
Cable Length	1.0m or 2.0m
Cable Type	PFA Flat Pair with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC