

## Standard-Thermoelement-Steckverbinder Inline-Buchse IS-K-F Typ K IEC



**Standard in-line Socket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard-Thermoelement-Steckverbinder Inline-Buchse IS-K-F Typ K IEC**

Nimmt Standard-Rundsteckerstifte auf, geeignet für industrielle Anwendungen

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-1402-001
<b>General Description</b>	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Standard
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS