




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## JIS Typ K Thermoelement mit Magnetleiste, PTFE-isoliert mit Miniaturstecker - 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-0554-001
<b>General Description</b>	Ideal for ferrous metal surface measurement
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Magnet</b>	50 x 25 x 1mm magnetic strip
<b>Cable Length</b>	2M
<b>Cable Type</b>	PFA insulated wires, twin twisted
<b>Termination Type</b>	Miniature flat pin plug
<b>Max. Temperature</b>	100°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Colour</b>	Positive: Red / Negative: White