




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Typ K Thermoelement mit Magnetleiste, PTFE-isoliert mit Miniaturstecker - 2m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0554-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	2M
Cable Type	PFA insulated wires, twin twisted
Termination Type	Miniature flat pin plug
Max. Temperature	100°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)
Colour	Positive: Red / Negative: White