




JIS Silikonkautschuk Patch Thermoelement, PFA Twin Twisted mit Miniaturstecker oder blanken Schwänzen - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

General Description	Resistant to many oils & chemicals
Sensor Type	Thermocouple junction encapsulated in a silicone rubber patch
Thermocouple Type	K
Cable Length	1m, 2m, 5m
Cable Type	PFA insulated twin-twist wire 1/0.2mm
Termination Type	Miniature flat pin plug or Bare tails
Max. Temperature	+150°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)
Colour	Positive: Red / Negative: White
Dimensions	Length 40mm x Width 13mm x Height 5mm Rubber patch