



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Feindraht-Scheiben-Thermoelement, PFA mit Bare-Tails-Anschluss - Typ K JIS - 1m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

Product Code	XE-0564-001
General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc, 6.35mm diameter x 0.3mm thick
Thermocouple Type	K JIS
Thermocouple Junction	Grounded measuring junction
Cable Length	1m
Cable Type	1/0.2mm PFA Teflon® insulated twisted wire
Termination Type	Bare Tails
Max. Temperature	200°C
Min. Temperature	0°C
Standards Met	JIS (Japanese Industrial Standards)
Colour	Positive: Red / Negative: White
Accessories	Supplied with high temperature adhesive kapton strips (13x18mm x 5) for quick and easy attachment to surface/process