




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Typ K Magnetstreifen-Thermoelement, PFA-isoliert mit Miniaturstecker



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

| | |
|------------------------------|---|
| General Description | Ideal for ferrous metal surface measurement |
| Thermocouple Junction | Grounded measuring junction |
| Magnet | 50 x 25 x 1mm magnetic strip |
| Cable Length | 1 or 2 metre |
| Cable Type | PFA insulated wires, twin twisted |
| Termination Type | Miniature Plug |
| Standards Met | ANSI MC96.1 |
| Colour | Positive: Red / Negative: Yellow |