





LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Typ K Magnetstreifen-Thermoelement, PFA-isoliert mit Miniaturstecker



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

General Description	Ideal for ferrous metal surface measurement
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	1 or 2 metre
Cable Type	PFA insulated wires, twin twisted
Termination Type	Miniature Plug
Standards Met	ANSI MC96.1
Colour	Positive: Red / Negative: Yellow