


ANSI Typ K Magnetstreifen-Thermoelement, PFA-isoliert mit Miniaturstecker



[Click here to view product](#)

Product Code	XE-0579
General Description	Ideal for ferrous metal surface measurement
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	1 or 2 metre
Cable Type	PFA insulated wires, twin twisted
Termination Type	Miniature Plug
Standards Met	ANSI MC96.1
Colour	Positive: Red / Negative: Yellow