



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **Bolzen-Thermoelement, Glasfaser-Edelstahl ber geflochtenem Kabel - JIS Typ K - M8x1mm Gewindesteigung - 13mm Bolzengre**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bolzensensoren werden in der Regel in die Oberfläche eingeschraubt, was in der Regel die Reaktionsgeschwindigkeit und die Empfindlichkeit der Messung erhöht. Geeignet für den Einsatz bis +250°C.

### Leistungsbeschreibung

## Specifications

<b>Product Code</b>	XE-0590-001
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Bolts have grounded junctions
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel over braided cable
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Bolt thread/pitch</b>	M8 x 1mm
<b>Colour</b>	Positive: Red / Negative: White
<b>Length</b>	13mm