



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Bolzen-Thermoelement, Glasfaser-Edelstahl ber geflochtenem Kabel - JIS Typ K - M12x1mm Gewindesteigung - 13mm Bolzengre



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bolzensensoren werden in der Regel in die Oberfläche eingeschraubt, was in der Regel die Reaktionsgeschwindigkeit und die Empfindlichkeit der Messung erhöht. Geeignet für den Einsatz bis +250°C.

Leistungsbeschreibung

Specifications

Product Code	XE-0594-001
Thermocouple Type	K JIS
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2m
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare Tails
Max. Temperature	250°C
Standards Met	JIS (Japanese Industrial Standards)
Bolt thread/pitch	M12 x 1mm
Colour	Positive: Red / Negative: White
Length	13mm