

**ANSI Typ K 1,9 kg
Zugknopf-Magnet-Thermoelement, PFA-isoliertes
Kabel mit Edelstahl-bergeflecht mit Miniatur- oder
Standardstecker**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

General Description	1.9kg Pull Button Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height
Cable Length	1.0m or 2.0m
Cable Type	Glassfibre insulated lead with stainless steel overbraid 7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	ANSI