

ANSI Typ K 1,9 kg Zugknopf-Magnet-Thermoelement, PFA-isoliertes Kabel mit Edelstahl-bergeflecht mit Miniatur- oder Standardstecker



[Click here to view product](#)

Product Code	XE-0620
General Description	1.9kg Pull Button Magnet Thermocouple
Thermocouple Type	K
Magnet	19.1mm Diameter / 12.7mm height
Cable Length	1.0m or 2.0m
Cable Type	Glassfibre insulated lead with stainless steel overbraid 7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	250Å°C
Min. Temperature	-50Å°C
Standards Met	ANSI