



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Gefertigtes geerdetes Thermoelement aus  
Edelstahlrohr, berflochtenes  
Glasfaser-Edelstahlkabel - Typ K JIS - 3,18 mm  
Durchmesser x 25 mm Lnge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

<b>Product Code</b>	XE-0644-001
<b>Sensor Type</b>	316 stainless steel tube in various sizes
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip for fast response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel overbraided cable
<b>Cable Termination</b>	Bare Wire Tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Dimensions</b>	3.18mm Diameter x 25mm length