



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Gefertigtes geerdetes Thermoelement aus Edelstahlrohr, beflochtenes Glasfaser-Edelstahlkabel - Typ K JIS - 4,0 mm Durchmesser x 12 mm Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-0646-001
Sensor Type	316 stainless steel tube in various sizes
Thermocouple Type	K JIS
Thermocouple Junction	Grounded at tip for fast response
Cable Length	2m
Cable Type	Glassfibre stainless steel overbraided cable
Cable Termination	Bare Wire Tails
Max. Temperature	350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	4.0mm Diameter x 13mm Length