



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## **ANSI Typ K 9kg Zugkraft (Hufeisen) Magnet-Thermoelement, PFA-isoliertes Kabel mit Edelstahl-bergeflecht mit Miniatur- oder Standardstecker**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded at tip
<b>Magnet</b>	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
<b>Cable Termination</b>	1.5 or 3 metres PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
<b>Termination Type</b>	Fitted Mini or Standard Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	ANSI MC96.1