


**JIS Typ K 9 kg Pull Power (Hufeisen)  
Magnet-Thermoelement, PFA-isoliertes Kabel mit  
Edelstahl-bergeflecht mit Abschluss im  
Standardstecker - 1,5 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Leistungsbeschreibung

**Specifications**

<b>Product Code</b>	XE-0680-001
<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Sensor Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip
<b>Magnet</b>	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
<b>Cable Length</b>	1.5m
<b>Cable Type</b>	PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
<b>Termination Type</b>	Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)