

**JIS Typ K 9 kg Zugkraft (Hufeisen)
Magnet-Thermoelement, PFA-isoliertes Kabel mit
Edelstahl-bergeflecht mit Miniatur- oder
Standardstecker**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded at Tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Termination	1.5 or 3 metres PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
Termination Type	Fitted miniature or standard plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)