


Standard-Rechteck-Thermoelement-Steckverbinder-Blendenbuchse IS-K-FF Typ K IEC



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Rechteck-Thermoelement-Steckverbinder-Blendenbuchse IS-K-FF Typ K IEC

Nylon-Clip zur schnellen und einfachen Befestigung an der Platte. Alle Blendenmontagebuchsen sind mittels ungleich groer Stifte polarisiert

Leistungsbeschreibung**Specifications**

Product Code	XE-1407-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	K IEC
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Green
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Standard Connector