




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Blatt-Thermoelement, glasfaserisoliertes Kabel mit Edelstahl-bergeflecht - Typ K JIS mit Standard-Steckeranschluss - 2m Lnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

Product Code	XE-0699-001
General Description	Thermocouple mounted on a stainless steel shim
Thermocouple Type	K JIS
Thermocouple Junction	Junction grounded at tip for fast response
Cable Length	2M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	25x13mm shim (0.2mm thick)