


Standard rechteckiger Thermoelement-Steckverbinder Faszienbuchse IS-N-FF Typ N IEC



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Rechteck-Thermoelement-Steckverbinder-Blendenbuchse IS-N-FF Typ N IEC
 Nylonclip zur schnellen und einfachen Befestigung an der Platte. Alle Blendenmontagebuchsen sind mittels ungleich groer Pins polarisiert
Spezifikationen

Specifications

Product Code	XE-1557-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	N IEC
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Pink
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS